AGA Sanitätsartikel

AGA-OP-POWER-MAT Operation Tables

The operation tables AGA-OP-POWER-MAT stand out not only due to their functionality but also to their elegant anodized lifting column. The hydraulic height adjustment is done with a centrally located foot pump.

Accessories

14-19.

• Paper roll holder (PU-2049).

· Central locking castor system

• Retrofitting is not possible.

• Further accessories online

and on catalogue pages

(LR/Z-50, Ø 50 mm).

Features

- The lifting column is made of anodized light alloy.
- All upholstered parts are manually adjustable. Through this type of construction, which is independent of a power source, the AGA-OP-POWER-MAT is available at a very competitive price.
- The removable upholstered segments (thickness 65 mm) are made of soft polyurethane foam and covered with a black, electrically conductive synthetic leather (cotton/polyester mix and 100% PVC compound coating).
- The headrest is continuously adjustable (positive and negative) and secured with clamping rod locking mechanism.
- The back, seat and leg section adjustments are assisted by gas pressure springs. The leg rests can be individually folded down by 90° .
- The standard side rails on the outside frame of the backrest, seat and leg sections allow the use of various accessories such as anesthesia pole, instrument table, side and shoulder supports, arm supports etc.



4 MP-1050/5Y, LR/Z-50

AGA Sanitätsartikel GmbH Löhner Strasse 198b, D-32584 Löhne, Germany Tel: +49 (0)5731 82283 · info@agasan.com



1 MP-1050/5Y, LR/Z-50

Adjustment Ranges

- Headrest +20°/-30°
- Backrest +45°/-20°
- Trendelenburg 20°
- Reverse Trendelenburg 8°
- Seat section +40°/-10°
- Leg sections 0°/-90°
- Height 785 1085 mm

More info online and in the AGA product catalogue on page 8.



L 1900 mm 295 575 360 520 W 500 mm 500 mm Max. load 150 kg Prod. weight 110 kg

2 MP-1050/5Y, NB-30560, BKO-2522, LR/Z-50

Picture	Article Number	Width (mm)	Frame		Price
Picture			Stainless Steel	Colour	(EUR)
in 1 - 4	MP-1050/5Y	500		х	
similar 1 - 4	MP-3050/5Y	500	х		

agasan.com



AGA-OP-POWER-MAT Equipment Options

The operation tables AGA-OP-POWER-MAT have been designed for various fields of applications and they may be equipped with different components accordingly. Below is a description of the available equipment options. Thanks to the mechanical adjustments, the AGA-OP-POWER-MAT remains independent of a power source.



5 MP-1050/5YGYN, MP-RF, VA-30/MP, LR/Z-50, PU-2049

Options

- Plastic-coated steel frame and plastic base cover (Article No. MP-1050/5Y).
- Glass bead blasted stainless steel frame and stainless steel base cover (Article No. MP-3050/5Y).
- Compartment for X-ray tape (Article No. MP-RF), approx. 435 x 305 mm, under the seat section. Retrofitting is not possible.
- Stainless steel tray (Article No. VA-30/MP), approx. 250 x 190 x 45 mm. Retrofitting is not possible.
- Upholstered foam pads (Article No. MP-PAL), thickness 40 mm. The upper surfaces are seamless. Hooks and loops fasteners are located on the undersides to attach the pads to the 15 mm thick radiolucent support plates.
- Removable, split and vertically/laterally adjustable leg plates (Article No. MP-BTSA) instead of the standard leg sections which can only be dropped down 90° (Vertical adjustment: 0°/-90°, lateral adjustment/spread: 0°/90°).
- The AGA-OP-POWER-MAT is also available with an adjustable headrest by gas pressure spring (Article No. KVG), instead of headrest adjustment with clamping rod locking mechanism.

Set combinations:

- Article No. MP-1050/5YGYN includes the operation table MP-1050/5Y with gynecological cutout, MP-PAL and MP-BTSA.
- Article No. MP-3050/5YGYN is equipped identically to MP-1050/5YGYN, except for the glass bead blasted stainless steel frame and stainless steel base cover.





6 MP-1050/5YGYN, MP-RF, VA-30/MP, LR/Z-50, PU-2049 7 MP-1050/5YGYN, MP-RF, VA-30/MP, BHG-600, BKG-600, RL-2060, LR/Z-50, PU-2049

Accessories

- Göpel leg rests (Article No. BHG-600) and Göpel mounting clamps (Article No. BKG-600)
- Ring armrests (Article No. RL-2060)

Further accessories can be found online and on catalogue pages 14-19

Picture	Article Number	Width (mm)	Frame		Price
			Stainless Steel	Colour	(EUR)
in 5 - 7	MP-1050/5YGYN	500		х	
similar 5 - 7	MP-3050/5YGYN	500	х		

AGA Sanitätsartikel GmbH Löhner Strasse 198b, D-32584 Löhne, Germany Tel: +49 (0)5731 82283 · info@agasan.com

agasan.com