



TRUSTED + QUALITY

The ergonomic solution for improved control

Injector Body

- + Crystal clear material for improved visibility during the delivery process
- + Designed for both mono and bi-manual injection procedure
- + Body shape and soft-touch texture for safe and comfortable handling

Spring

- + Allows push back for safety
- + Enhanced feeling for controlled delivery

Cartridge


- + Rounded corner for safe cartridge loading
- + Available in two versions: GMS and Coated
- + Cartridge sizes from 1.8 mm to 2.6 mm

Plunger

- + Designed with a safety stopper for a controlled start
- + Engineered design and robust material for secure alignment

Soft Touch

- + Innovative anti-slip texture for safety (when using wet gloves)
- + Comfort of use





TRUSTED⁺ QUALITY

The ergonomic solution for improved control

For a controlled delivery

- + Injector body design for better handling and made of crystal clear material with large fenestrations for improved visibility.
- + Robust plunger designed with an advanced material for perfect alignment and a stopper for safe stop/start.
- + Soft touch and anti-slip texture for safe and comfortable handling.

Available in two different cartridge versions

ErgoTouch and ErgoTouch Plus have a safe closure mechanism and are available in a wide range of cartridge sizes from 1.8 mm to 2.6 mm to suit the chosen surgical technique.



- + GMS gliding agent
- + Controlled residual levels
- + Established technology with standard PP



- + Unique coating treatment technology with special PP
- + No residues with no additive transfer