

LOJER[®] 100
Since 1919

Future disinfectant recommendation



Standard colour



Standard colour



Standard colour



Future disinfectant recommendation

Future	Disinfectant	Manufacturer
XX	Antibac	Antibac AS
XX	Meliseptol foam pure Hexaquart	B. Braun Melsungen AG
XX	Perasafe	Berner Pro
XX	Virkon	Berner Pro
XX	Kohrsolin	BODE Chemie GmbH
XX	Desomed rapid	Desomed Dr. Trippen GmbH
XX	Cleanisept Wipes	Dr. Schumacher GmbH
XX	FD360	Dürr Dental AG
XX	FD360	Dürr Dental AG
XX	Brial	Ecolab
XX	Klercide 60/40	Ecolab
XX	Incidin Oxywipe	Ecolab
XX	Easydes	KiiltoClean Oy
XX	Klorilli	KiiltoClean Oy
XX	Surfanios	Laboratoires Anios
XX	Surfa'safe	Laboratoires Anios
X	Aldasan	Lysoform Dr. Hans Rosemann GmbH
X	Antifect plus	Schülke & Mayr GmbH
-	Buraton rapid	Schülke & Mayr GmbH
XX	Microzid	Schülke & Mayr GmbH
XX	Wet wipe	Wet Wipe A/S

xx = recommended

x = recommended to use diluted, 0.50% dilution

- = not recommended